

## Goat anti-Rabbit IgG (H+L) Secondary Antibody, DyLight<sup>™</sup> 488

Product Details	
Size	100ul
Species Reactivity	Rabbit
Host/Isotype	Goat / IgG
Class	Polyclonal
Туре	Secondary Antibody
Conjugate	DyLight™ 488
Immunogen	Purified Rabbit IgG, whole molecule
Form	Liquid
Concentration	1 mg/mL
Purification	Antigen affinity chromatography
Storage buffer	PBS, pH 7.2, with 1% BSA
Contains	0.02% sodium azide
Storage conditions	4° C

Applications	Tested Dilution	Publications
Western Blot (WB)	1:1000-1:10000	-
Immunohistochemistry (IHC)	1:50-1:2000	-
Immunocytochemistry (ICC/IF)	1:50-1:2000	-
Flow Cytometry (Flow)	1:25 - 1:100	-
Immunoprecipitation (IP)	Assay-dependent	-

DyLight 488 Amax= 493 nm; Emax= 518 nm. Mole Dye/Mole Protein Ratio is lot-dependent.

For Research Use Only. Not for use in diagnostic procedures. Not for resale withort express authorization, Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of a set. as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications is made. This warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, the warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTES, EXPRESS OR IMPLED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLED WARRANTES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR UNROSE, CREVEND FOR NON INFRINCEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFONING PRODUCTS DURING THE WARRANTES OR REPUND FOR REPUND FOR REPUND FOR REPUND FOR REPUND FOR THE NON-CONFORMING PRODUCTS), STELE ESS SUCE OFTICE, SS SUCE OFTICE, SS SUCE OFTICE, SS SUCE OFTICE, SS SUCE OFTICE, ST SUCE OFTICE, SS SUCE OFTICE, SS SUCE OFTICE, ST SUCE OFTICE, SS SUCE OFTICE, ST SUCE OFTICE, SS SUCE OFTICE, SS SUCE OFTICE, ST SUCE OFTICE, SS SUCE OFTICE,

Note: the product information may be optimized and upgraded. Please refer to the actual label information.